
SkinElly

INTENSE SUB-Q

Hyaluronic acid gel

INSTRUCTION FOR MEDICAL USE OF MEDICAL PRODUCT

Medical Product Name:

Intradermal Implant «SkinElly» with Sodium Hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml composed of:

- 2 ml syringe with an intradermal implant with sodium hyaluronate - 1 pc.
- Needles - 0,5x16 mm (25 G)-2 pcs.
- Self-adhesive labels (sticker) - 3 pcs.
- Instruction for use (package leaflet) – 1 pc.

MANUFACTURER/DEVELOPER OF THE MEDICAL PRODUCT:

ACTIVE MEDICAL PRODUCTS LIMITED, Suit 1,
5 Percy street Fitzrovia, London W1T 1DG (UK)

CORRECT USE OF THE MEDICAL PRODUCT WITH INDICATION OF THE POTENTIAL CONSUMER:

Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml is indicated for use in cosmetology, plastic and aesthetic surgery for the correction of involuntal skin changes. Is to be injected into the dermis (depending on the approach).

Qualified prospect of a medical product - Intradermal implant with sodium hyaluronate called "SkinElly", patients with paracmastic skin changes.

Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml is an injection, transparent, colorless, sterile, pyrogen-free, fully biosynthetic, 100% bioabsorbable in human body tissues, placed in a sterile glass syringe with a volume of 2.0 ml pre-filled with a colorless, transparent and viscoelastic substance.

INGREDIENTS:

Hyaluronate sodium 24 mg/ml (2,4 %)

Phosphate buffer solution

Diglycidyl ether 1,4-butanediol <2 ppm

INDICATIONS FIELD OF THE MEDICAL PRODUCT:

Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml is used in plastic and aesthetic surgery, aesthetic medicine and cosmetology to correct involuntal skin changes. Is to be injected in deep dermis supraperiostal to make a volumetric modeling of the middle and lower third of the face, correction of gravitational presentations of degrees 2 or 3. It is used

to correct expressed deficiency of tissues, to recreate the volume of the cheek and cheekbone area, clear contour lines of the facial contours, chin and jaw angles augmentation. Ideal for patients with deformation and tired type of aging and expressed gravitational ptosis, with expressed volume losses.

CONTRAINDICATIONS TO THE MEDICAL PRODUCT USE:

General:

- Medical products only for intradermal use;
- Allergic reaction to any components that make up the medical product;
- Pregnancy and lactation;
- In combination with laser procedures or peeling;
- Blood vessels administration;
- Systemic treatment with retinoids;
- Systemic anticoagulant and antithrombotic therapy;
- Systemic immunosuppressive therapy;
- Inflammatory skin infiltration;
- Progressive chronic dermatosis (Koebner Phenomenon);
- Infection caused by the herpes virus in the acute stage;
- Paraneoplastic process at the administration site;
- Severe somatic disease at the decompensation stage;
- Acute infectious diseases, feverish conditions;
- Acute exacerbation of chronic diseases, accompanied by decompensation processes of the II-III stage;
- Autoimmune diseases;
- long-term treatment at the dentist (dental prosthetics) during the procedure;
- Under the age of 18 (without the consent of the parents or their substitute);
- Hypersensitivity to hyaluronic acid drugs in past medical history;
- Bleeding disorder: hemophilia;
- Porphyria.

Absolute:

- Tendency to form hypertrophic and keloid scars or any other impaired healing;
- Permanent fillers in the administration area.

Relative:

- Inflammation and infectious processes in the intended administration area (acne, herpes, eczema, etc.);
- Violations of the epidermis integrity in the administration zone (for example, post-pilling period, etc.).

PRECAUTIONS:

There were no studies on the use of a medical product depending on the form of execution of the Intradermal implant with «SkinElly» sodium hyaluronate in patients

with already implanted medical products based on hyaluronic acid. It is necessary to determine the "benefit – risk" ratio of the proposed procedure before using medical product.

- Check the expiration date and integrity of the medical product packaging before use;
- Use only for one patient;
- Do not re-sterilize;
- Comply with sanitary and epidemiological requirements during injection procedures as there is a risk of infection;
- Do not expose to intense heat (sunshine or tanning beds) or significant cooling of the administration site;
- Dispose of the syringe with implant residues after use.

Note: Removing or replacing an **intra dermal implant** is not possible. The use of various hardware techniques (for example, d'Arsonval, electroporation, microcurrent therapy, vaporization) or the administration of the enzyme hyaluronidase or Longidase is possible for removal or accelerated absorption at the discretion of the practitioner.

INTERACTIONS WITH OTHER MEDICAL PRODUCTS AND MEDICINAL DRUGS

The use of a medical product Intra dermal implant «SkinElly» with sodium hyaluronate in combination with other medical products and medicinal drugs has not been studied.

Note 1: Removing or replacing an Intra dermal implant «SkinElly» with sodium hyaluronate is not possible. The use of various hardware techniques (for example, d'Arsonval, electroporation, microcurrent therapy, vaporization) or the administration of the enzyme hyaluronidase or longidase is possible for removal or accelerated absorption at the discretion of the practitioner.

Note 2: The effect of the Intra dermal implant «SkinElly» with sodium hyaluronate on vision and reaction time has not been analyzed and therefore there are no special instructions for drivers and/or persons operating complex technological devices.

The doctor should not:

- Use together with other drugs, both in one syringe, and in one procedure;
- Administer to patients with known bleeding disorders or patients receiving medicinal drugs that increase blood clotting time (vitamin K antagonists, heparin, or immunosuppressant drugs), as they may increase the risk of bleeding at the site of administration;
- With unsatisfactory experience of this type of injection;
- Use before or after laser treatment or chemical peeling, as this may cause inflammation of the injected area.

The patient should not:

- Subject the treated area to heat (sunbathing, sauna, steam bath, etc.) or extreme cold until any signs of local inflammation disappear.

There is a risk of infection characteristic of all subcutaneous injections. Therefore, standard precautions are used for this type of procedure when administering the product.

You must ensure that the individual sterile packaging is not damaged, the expiration date of the syringe and needles has not expired before use.

Do not use the product if the contents of the syringe look cloudy.

It is forbidden to administer into blood vessels, bones, ligaments, tendons and muscles and into birthmarks.

The remains of the medical product must be disposed of!**It is forbidden to use medical product:**

- With sterile packaging integrity breach. Do not re-sterilize!
- If the syringe cap of the Luer-Lok cap is gapping;
- Expired;
- With color changes or transparency of the syringe contents.

**EXPECTED AND PREDICTABLE
ADVERSE REACTIONS****Transient reactions:**

As a manifestation of transient reactions, the appearance of short-term bleeding (1-2 minutes) and the appearance of bruising, pain and itching at the moment of injection for a short time (10-15 minutes) are possible. Perhaps the appearance of erythema, as a result of local irritating or compensatory action, which, in turn, passes quite quickly (no more than 2 days). The manifestation of the injection site mass (for no more than 4-5 days), discoloration of the skin and an increase in sensitivity at the injection site.

Preservation of any reactions over a specified time and (or) the addition of inflammatory symptoms (temperature, infiltration, fluctuation, puffiness and swelling, granuloma, etc.) indicate complications and require that the patient be treated by a specialist with appropriate profiles using pharmaceuticals drugs.

BIODEGRADATION

Biodegradation of the Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml depends on the design (concentration) and the individual characteristics of the patient:

BIODEGRADATION OPTIONS

Hydrolytic degradation: The destruction of the material in the aquatic environment under the influence of:

- acids and alkalis (non-enzymatic hydrolysis);
- enzymes (enzymatic hydrolysis);
- reactive oxygen intermediates and free radicals (oxidative destruction, which is significantly accelerated in the presence of metal ions as a catalyst).

Cell degradation: Biomaterial destruction by macrophages and mast cells.

Bacterial degradation: The presence of bacteria in the area of implant insertion can accelerate biodegradation due to enzymatic and oxidative degradation.

Mechanical degradation: Fragmentation and/or change of the implant shape under the action of static or dynamic loads can accelerate the destruction of the material (active facial expressions, massage, strokes, etc.).

ELIMINATION PATHWAYS:

Sodium hyaluronate is eliminated from the body depending on the individual characteristics of the patient; in the human body sodium hyaluronate is laminated by hyaluronidase enzymes. The decomposition products of sodium hyaluronate are oligosaccharides and extremely low molecular hyaluronates, which are naturally eliminated from the human body through the liver and kidneys.

TIME OF ONSET OF THERAPEUTIC/ COSMETOLOGICAL EFFECT

The time of onset of the therapeutic/cosmetological effect depends on the individual characteristics of the patient.

The time of therapeutic effect preservation

The time of therapeutic effect preservation depends on the individual characteristics of the patient and in the human body.

Information on the procedure of working with a medical product (procedure and/or usage rules).

Preparation for Use

The procedure is carried out in accordance with the standard medical practice and local regulations. The product should be used by doctors who have received special training in injection techniques for filling skin folds for a successful procedure.

Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml to be used as it is delivered.

Carefully inspect all parts for damage.

- Before injection, carefully inspect the SkinElly implant visually. You may need to perform several additional procedures to achieve optimal results. In order to enable effective evaluation of the implantation results, a break should be made between the courses of treatment for at least seven days.
- Before starting treatment, patients should be informed about the indications for the product use, contraindications, incompatibility and potential side effects. Before the procedure, you should read the full history of the patient and carefully examine the area to be treated.
- You should also evaluate the patient's ability to cope with pain and, if necessary, apply the most appropriate form of anesthesia. In the case of administration of an anesthetic substance, ice should be applied to the treated area to reduce swelling and local stretching.
- Before starting the procedure, thoroughly clean the area to be treated with soap and water and disinfect with an alcohol swab.

Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml: is to be injected supraperiosteal in deep dermis.

Injection zones:

Face, neck, décolleté, hands, earlobes, limb skin, trunk skin.

Injection technique:

- Micropabules.
- Microdroplet to the area of "kiset" wrinkles, "proud" wrinkles.
- Microbolus in the red border of lips, earlobes.
- Linear-retrograde with separation of superficial linear wrinkles.

Immediately stop the procedure if you suspect a puncture of the vessel.

If significant resistance is felt during the gel injection, the needle should be moved deeper to ensure easier insertion of the implant. If resistance remains, try inserting the implant from another entry point. In case of continuous resistance, it may be necessary to replace the injection needle and/or syringe.

Any excessive pressure on the syringe should not be applied at any time during the procedure.

If blanching occurs immediately after the gel injection, discontinue the injections and massage the treated area until the skin returns to its normal color. If the skin has not regained its normal color, the product administration process should not be resumed.

Superficial injections or the injection of large doses of the product may cause discoloration and nodulation.

The duration of the correction of the Intradermal implant with sodium hyaluronate SkinElly depends on the choice of performance, the area of administration and the characteristics of the patient. During the procedure after the injection, a light massage is desirable for better redistribution in the tissues.

INSTRUCTIONS FOR THE PATIENT

- Patients should adhere to the following recommendations: avoid excessive physical exertion, as well as exposure to direct sunlight and tanning lamps or extreme weather conditions for 24 hours.
- Patients should be informed of the application of a wrap with ice or a cold compress to the treated area within 24 hours after the procedure.
- Patients are also encouraged to talk less, smile and laugh for one week after the procedure.
- Patients who underwent oral procedures should rinse mouth with saline every 3–4 hours, within one week after the procedure.
- If necessary, patients can use oral analgesics as recommended by their physician.
- If the inflammatory process lasts more than a week or other side effects appear, you should immediately consult a doctor.
- You should refrain from applying makeup for 12 hours after the injection, you should avoid prolonged contact with sunlight, ultraviolet light, gels, do not go to the sauna or the bath for 2 weeks after the injection.
- Do not consume alcoholic beverages, anti-inflammatory drugs, anticoagulants and drugs accelerating the metabolism before and after the procedure without consulting a doctor.

THE OPERATING CONDITIONS OF THE MEDICAL PRODUCT:

The Intradermal implant SkinElly with sodium hyaluronate can only be used by medical personnel with appropriate qualifications and familiar with injection techniques.

- Temperature from + 32 °C to +42 °C;
- Humidity from 10 to 80% (non-condensing);
- Atmospheric pressure from 500 to 1060 hPa.

Avoid direct sunlight, freezing, and condensation.

INCOMPATIBILITY

Since hyaluronic acid can cause the release of a disinfectant tertiary (medium) ammonium, for example, benzalkonium chloride or chlorhexidine, the drug should not be near these substances or come into contact with the devices with which these substances were used.

Upkeep and periodic maintenance, including cleaning and disinfection

In the production process of SkinElly, hyaluronic acid undergoes several steps of purification from impurities of residual proteins and endotoxins peculiar to bacterial cell.

Due to the low concentration of residual ingredients of bacterial origin, SkinElly does not exhibit allergenic properties and allergy test before the procedure is not required.

STORAGE CONDITIONS AND TRANSPORTATION STORAGE AND ACTION AFTER USE

After completion of treatment, the syringe must be disposed of in accordance with established rules. You must dispose of the syringe with the remnants of the gel and the needle, as they are potentially infected.

STORAGE CONDITIONS

For Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml

- Temperature from + 2 °C to +25 °C;
- Relative humidity 35-85%;
- Atmospheric pressure 500 - 1060 hPa.
- Non-condensing, away from sunlight.

Ensure packaging integrity before use.

TRANSPORTATION CONDITIONS

Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q Version 24 mg/ml the temperature is + 2 to +25 °C

DISPOSAL OF A MEDICAL PRODUCT

Disposal of a medical product Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly INTENSE SUB-Q version 24 mg/ml in the territory of the Russian Federation is carried out according to SanPiN 2.1.3684-21 belongs to Class B (epidemiologically hazardous waste) and is a potentially infectious waste.

The used medical product, including needles for injection, not fully used syringes are disposed of according to the specified hazard class of medical waste according to SanPiN 2.1.3684-21. Expired medical products, medical products in damaged packaging should be disposed as class A of medical waste according to SanPiN 2.1.3684-21.

THE SHELF LIFE OF THE MEDICAL PRODUCT

The shelf-life is 2 years from the moment of sterilization while meeting the conditions of transportation and storage


The opened medical product should be used immediately and disposed of immediately after use.

OFFICIAL REPRESENTATIVE CONCERNING THE MEDICAL PRODUCTS IN THE RUSSIAN FEDERATION

ACTIVE MEDICAL PRODUCTS Limited Liability Company

Russia, 123100, Moscow, 16 Shmitovsky Passage, bld. 2,
Room 1, Office 13, Floor 2

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